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
<使用说明书>

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
# Folding 7-color Spectrometer

<Instructions>

In order to use this product correctly,
please read the manual carefully and keep it safe after reading.

# Safety precautions

**Before use, it is necessary to read the operating manual and use it correctly.**

The user manual not only provides instructions for the use of this device, but also includes safety precautions, storage, and product warranty information. For safety reasons, please confirm all the following safety information and start-up instructions before using this device. We are not responsible for any personal injury or material loss caused by non-compliance with the above clause.

## General considerations

Before using the equipment, it is necessary to inspect the body and LED module to be used.
This device should be kept away from high temperatures, high humidity, and direct sunlight, and should be placed indoors with good ventilation.
- Recommended temperature during startup: 0-26 ℃
- Recommended temperature for storage: 0-35 ℃
- Relative humidity: below 80%

Please set it indoors with good ventilation and avoiding direct sunlight.
Please install while maintaining a level surface in a relatively flat area.
This device may also cause significant damage under weak impacts, so do not apply external impacts. Be sure to move it slowly during handling.

## Notes on the hazards of power supply

It is prohibited to use the device simultaneously with other heat transfer products during startup.
Please be careful not to let water enter the equipment, as the high voltage current formed by moisture can cause fire or electric shock.
Please do not use damaged power plugs or loose sockets.
When not in use for a long time, unplug the power plug and clean it before storing it. Otherwise, it may cause electric shock, leakage, and fire.
Please unplug the power and contact the service center in a timely manner when there is water ingress, abnormal sound, burnt smell, smoke or other abnormal equipment in the product.


## Taboos


Even if it does not fall under the following precautions, please do not use it if instructed by the operator.
-Patients with heart disease, mental illness, and hypertension
-When suffering from skin diseases (infectious)
-When there is a disease similar to skin inflammation

## Regulatory matters

Users must be familiar with the user manual before using the device.
When using unapproved decorations, it may trigger device malfunctions, and any resulting damage will not be provided free of charge. Therefore, please use genuine accessories.
Please do not disassemble or modify maintenance equipment without authorization. Any issues caused by unauthorized disassembly/repair will not be compensated.

# Product composition


Adapter

Patch

## Product Settings

1.Before installation, it is necessary to consider the position and space of the equipment and maintain a certain distance from the heating equipment and other medical devices.
2. Avoid direct sunlight and high humidity, and install in a well ventilated area that can maintain suitable temperature and humidity.


Unfold one side panel  Unfold the other side panel

## Install

1) Open the packaging, gently remove the instrument body, unfold both sides of the instrument, and place it on a stable surface.
2) Remove the adapter power supply, insert the DC terminal into the instrument's DC connector, and ensure that the terminal connection is secure.
3) Insert the AC plug of the power adapter into the mains socket. At this time, the power is turned on and the instrument enters standby mode.
4) When not turned on, the time display text on the local panel is invisible.

# Instructions for use


## ON/OFF

After powering on, enter "standby mode"; short press this button to enter "power on mode", the indicator light will light up and display the default time of 30 minutes;
Press again to exit the power on mode and return to standby mode.
Note: When turning on the device, short press once to turn it off, and short press again to restart.


## Spectral switching

After turning on, short press this button to enter spectral mode, which has 8 modes to switch between. When the display shows 00, it is in standby mode.


1) Red light mode, display 01 on the monitor
2) Blue light mode, display 02 on the monitor
3) Green light mode, display 03 on the monitor
4) Yellow light mode, display 04 on the monitor
5) Orange light mode, display 05 on the monitor
6) Purple light mode, display 06 on the monitor
7) Cyan light mode, display 07 on the monitor
8) Loop light mode, display 08 on the monitor


1. Red light (wavelengths 647-760nm)


2. Blue light (wavelength 430-455nm)


3. Green light (wavelengths 492-550nm)


4. Yellow light (wavelength 550-580nm)


5. Orange light (wavelength 588-647nm)

6. Purple light (wavelength 360-430nm)

7. Cyan light (wavelength 480-490nm)


8. Circular light

## Time adjustment

The default startup time is 30 minutes. Press this button to increase the time in increments of 5 minutes, which can be continuously increased up to a maximum of 60 minutes; If you continue to press this button, it will jump to 00 minutes of working time. The time adjustment range is 00-60 minutes. If there is no button operation, the instrument will automatically exit the set state and start working at the current set working time after 3 seconds.
Note: When the time is 00, the instrument stops working.

## Spray switch

The default power on state is off. Press this button to start the spray function.
There are three gears to select (01 is the lowest gear, 03 is the highest gear).
Note: When the display is 00, turn off the spray function.

## Ultrasonic fine spray

Use ultrasound to generate nano unit spray, tiny water particles penetrate deep into the skin, and the absorption speed of water supplement is faster; It is recommended to use diluted products or purified water, and it is forbidden to use other high concentration make-up water or essence oil in combination; It is better to use it with facial mask.

## Heating switch

The default power on state is off. Press this button to activate the hot compress mode, and the display will show 01 with a default temperature of 45 ℃.
Note: When the display is 00, turn off the hot compress mode.

## Hot compress mode

Hot compress is a form of physical therapy, and it plays an important role in soft tissue injury diseases. By warming the skin deeply, opening pores makes it easier to absorb nutritional products. Two to three times a day, for fifteen to twenty minutes each time.

Note:
When using hot compress, the local skin should be checked for changes at all times. If there is any discomfort, it should be stopped immediately.

# Precautions during use

1) It is prohibited to directly look at any type of light emitted by this device with your eyes. When shining on the face, it is recommended to wear sunglasses or a blackout eye mask.

2) Patients with allergic skin should be cautious in conducting small-scale tests in concealed skin areas, and only use them on a large scale if there is no allergy.

3) If any discomfort occurs during the use of this product, it should be immediately stopped to recover. When using it again, the principle of short-term and frequent exposure should be followed. Wait until the body gradually adapts before using it normally.

4) This machine should be placed in a cool, dry, ventilated area that is not easily accessible to children. It is prohibited to place it in an environment with direct sunlight, humidity, high temperature, strong acid, strong alkali, strong oxidant, and strong dust.

5) Do not pour water onto the instrument; Do not use hands with water to operate this machine.

6) When using this device, please keep a certain distance from other electrical appliances, away from strong magnetic interference and flammable and explosive materials (gasoline, diesel, fireworks, etc.).

7) To clean this machine, please gently wipe it with a wrung out wet towel. Do not use strong acid or alkaline chemical reagents for cleaning. Soft brushes can be used to clean gaps or step positions.

8) Collisions, falls, throwing, and tapping of this device are prohibited, and any objects (including but not limited to mobile phones, bank cards/access cards, and other magnetic cards) are prohibited from being placed on the instrument.

9) It is prohibited to insert conductive objects such as screws and metal wires into the interior of this machine to avoid burning, short circuit, fire, and electric shock.

10) Minors and individuals without autonomous activity/perception/expression abilities should use this device with the care and assistance of a person with normal operation and perception abilities, or follow medical advice.

11) The use of unofficial accessories is prohibited, and modifications and self disassembly of this machine are prohibited. The above are not covered by the warranty.

# Product specifications

| | |
|---|---|
| Product name | Folding 7-color Spectrometer |
| Producer | China |
| Power supply | DC12V 3A 50HZ/60HZ |
| Power | About 36W |
| Product size | 480×220×337 (mm) |
| Product weight | 1.3kg ±10% |

# Maintenance

1. Clean the controller or body with a gentle towel.
2. When storing the product for a long time, cut off the power and avoid storing it in a place with high humidity.
3. When storing the product, keep the folding support folded up and store it.

Expand state

Folding state

Attention:
Users must cut off the power supply of the equipment and disconnect the power cord during maintenance and management.
Cleaning or storing the equipment while the power is turned on can lead to abnormal conditions.

## Matters to be confirmed before applying for faults

1.Before installation, it is necessary to consider the position and space of the equipment and maintain a certain distance from the heating equipment and other medical devices.
2. Avoid direct sunlight and high humidity, and install in a well ventilated area that can maintain suitable temperature and humidity.

| Problem points | Estimated reasons | Confirmation items |
|---|---|---|
| Abnormal power supply | Poor contact of adapter | Confirm power ON/OFF status<br>Confirm if the adapter is loose |
| Abnormal device startup | Poor LED status | Replace the corresponding LED panel when the LED is not lit(Contact purchasing agent/service center to exchange LED) |
| | Poor actuator status | Contact purchasing agent/service center |

# Maintenance

## Free warranty

1) Any malfunctions that occur during normal use within one year of purchasing the product are eligible for free service.
2) The free warranty only applies to the device body.

## Exclusions during the free warranty period

Without complying with the device startup steps, safety and management matters explicitly stated in this manual, you are unable to enjoy free services. Please advise
Pay close attention to the following:
1) When the adapter is pulled and causes damage to the device due to movement, it is excluded from the free warranty period.
2) The user's use of equipment or modification for other purposes is excluded from the free warranty period when disassembled.

## Paid warranty

1) After purchasing the equipment for one year, a paid warranty will be implemented, even if the equipment is damaged due to consumer negligence or other natural disasters (fire, salt flood, impact, earthquake, lightning, etc.) within the warranty period.
2) The paid guarantee items are as follows: damage to the device caused by consumer's unauthorized disassembly, modification, or careless use (water ingress, damage, etc.).

## Note

※This product is manufactured after a rigorous quality management and inspection process.

※When repairing for free, this certificate must be provided.

※If a malfunction occurs within the scope of normal consumer use, the manufacturer shall provide free repair services during the repair period.

※If the product is within the repair period but fails to use the original power supply, if the operator does not use it properly, if the instrument is dismantled, repaired, and the internal structure is modified on their own, or if the malfunction is caused by human damage or tapping, the manufacturer has the right not to provide repair and charge related component fees.

※The product guarantee letter will not be provided again, so it is essential to keep it properly.

| Model | | Warranty period | 1 year |
|---|---|---|---|
| Product name | | | |
| Purchasing date | | | |
| Corporate name | | | |
| Mailing address | | | |
| Full name | | | |

Made in China