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Soft rubber Spectrometer

Instructions

Thank you for purchasing this product! In order to use this product correctly, please read this manual before use;
In order to read the manual at any time when you need it, please keep it properly.

Safety precautions

Please read the safety precautions before use.
Please use the product correctly according to the instructions and suggestions.
The product shall be cleaned in time after use to avoid the accumulation of grease and its residues.

Danger	Contents related to personal injury caused by wrong operation
Be careful	Contents related to personal safety or property loss due to wrong operation

Danger

- Please use the connecting line supplied with the product, and connect the adapter plug to the household socket correctly, so as to avoid failure, fire or fire caused by heat.
- If the product is deformed, damaged or has peculiar smell, please do not use it again to avoid the risk of electric shock and fire.

Be careful

- The single illumination shall not exceed 30 minutes, and can be used twice a day with an interval of more than 6 hours.
- Do not drink, eat spicy and stimulating food, eat seafood and other high protein food within 24 hours after use, and avoid exposure to the sun (easy to tan).
- Please do not disassemble or modify the product by yourself, so as to avoid failure.
- If the product is damaged, in order to avoid danger, it must be repaired or replaced by professionals from the manufacturer, its maintenance department or similar departments.
- This product is not waterproof. Please do not wash it with water and take waterproof measures.

- Please do not store this product in humid environment.
- Do not use this product near inflammable and explosive materials, such as alcohol, gasoline, spray, etc., which may cause fire or explosion.
- Supervise children to make sure they don’t play with the product.
- Keep packaging materials such as plastic bags out of the reach of children.
- Before cleaning and maintenance or in case of failure in use, please shut down and unplug the adapter power plug.
- Avoid falling, hitting the adapter or squeezing the mask to avoid failure.
- Do not use this product for any purpose other than specified in this manual, so as to avoid failure or fire caused by heating.
- If you feel skin or body discomfort, please stop using immediately.

Do not use this product when:

— If you feel skin or body discomfort, please stop using immediately.

— The face has been filled or reshaped (hyaluronic acid and other implants), and the product cannot be used within one month.

— Patients with skin diseases, severe skin infection and other skin diseases.

— Patients with bleeding tendency and epilepsy.

— When you feel unwell or have a fever.

Introduce

This is a set of red light, blue light, green light, yellow light, blue light, purple light, laser light in one LED mask rejuvenation instrument. The product is composed of multiple LED lights, which can realize the light care of the face skin without dead corner coverage.

Component Description

Front

Opposite


Type–Cchargingport

USBoutputport

Display Description

Operating instructions

Controller

ON/OFF key：

- The default state is “OFF”, long press ON/OFF.
- After turning on, the default gear is 1st gear. Short press to adjust, 1–5 gears can be adjusted.

Color key：

- The default is “Red light mode”, “Green light mode” for the second time, “Blue light mode” for the third time, “Cyan light mode” for the fourth time, “Yellow light mode” for the fifth time, “Purple light mode” for the sixth time, “Laser light mode” for the seventh time, and “Cycle light mode”for the eighth time.

Time key：

- After turning on, the default time is 10 minutes.Short press to adjust, 1–90 minutes can be adjusted.

Cleaning and maintenance

- Please shut down and unplug the power plug before cleaning.
- Protect the appliance from moisture. Do not put the appliance in the dishwasher or microwave oven.
- The product is not waterproof and can not be washed with water. It can be cleaned with soft and wet cloth and wiped with soft and dry cloth.
- Do not use sharp or corrosive, volatile diluent cleaner.

Use operation

【Preparation before use】

- Clean the mask with a damp cotton pad.
- Cover the mask on the face and find a suitable and comfortable place to wear. The two sides of the mask pass through the ears to adjust the tightness of the magic belt to fix the mask.
- It is recommended to do photon cosmetology in the rest time as much as possible, find a comfortable sofa or bed, lie down quietly and comfortably, close your eyes and enjoy the beauty process brought by photons. Do not use phototherapy mask while reading books and magazines, watching TV, surfing the Internet, etc.

Use on a sofa or bed

You can't read or watch TV at the same time


Warm tips: users can choose different photon colors according to different needs, and do not use all three photons in a continuous time. According to the personal situation, different photons can be used continuously for beauty.

Specifications

Product Name: Soft rubber Spectrometer	
Voltage：5V3A	Rated power：20W
Product size：395x290x20mm	Net weight： 300g
Package size：465x320x35mm	

※Packing list: Host x1,Controllerx1,USB cablex1,Manualx1.

※Note: in view of the consistent policy of continuous improvement of products, the company reserves the right to change product specifications, and the physical products shall be subject to shipment.


Red light—LED

Wavelength:615–760nm

Blue light—LED

Wavelength:430–475nm



Green light—LED

Wavelength:492–550nm

Yellow light—LED

Wavelength:550–596nm


Cyan light–LED

Wavelength:455–492nm

Purple light—LED

Wavelength:360–430nm

White light—LED

After sales maintenance

- The appliance is only suitable for personal home care, not for commercial use.
- We provide a one–year warranty for our products, which starts from the date of purchase.
- During the warranty period, we will replace parts or the whole product as the case may be due to defects in materials or production process.
- Warranty claims must be submitted to the factory for after–sales service together with a certificate (purchase invoice) that can prove the date of purchase.

The company does not guarantee：

— Damage caused by improper use, such as using the wrong power supply, inserting the product into the improper plug, artificial damage, etc.

— The aging and wear of products caused by normal use have little effect on the use function of products.

— Product problems caused by improper operation or non professional maintenance or use of non company accessories for other purposes.

— Consumables, accessories and other consumables.

— Damage caused by irresistible factors (such as lightning, earthquake, fire, flood, etc.).

※ After sales service does not include damage caused by the above reasons.

Certificate

Inspection date:_____

Examination clerk: _____